

Readme

Group Project in Visualization of Medical Data

Jakob Troidl, 152501

Marian Stoschitzky, 1226251

Setup

1. Open the files **brain_visualization_2D.mlab** and **brain_visualization_3D.mlab** with MeVisLab (we used v3.2 for creation).
2. The data import in MeVisLab should work automatically, as we have used relative paths to the **_data_** sub-directory, where the imported files are stored. If that does not work, please update the import path in the modules *ImageLoad* and *vtkMINCImageReader* accordingly.

Visualizations

Double-click the “View2D” modules in the group of interest (*SumOfErrors*/*BoxPlotDeviation*) in the file of interest (**brain_visualization_2D.mlab/brain_visualization_3D.mlab**) for the corresponding visualization and dimensionality.

2D Visualization

You can change between coronal, axial and sagittal axis by (re-)wiring the corresponding output of the 3 *OrthoReformat3* modules to their succeeding modules (see figure 1).




Figure 1: Switch between *OrthoReformat3* outputs to switch between axes.